

Cebadores del Servicio de Secuenciación:

Código	Tipo
M13-20	M13 Seq Primer-20: 5'd[GTAAAACGACGGCCAGT]3' -17mer
M13R-24	M13 Rev Seq Primer-24: 5'd[AACAGCTATGACCATG]3' -16mer ¡Agotado!
M13R-24-17mer	M13 Rev Seq Primer-24: 5'd[CAGGAAACAGCTATGAC]3' -17mer
SP6	SP6 Universal Primer: 5'd[ATTTAGGTGACACTATAG]3' -18mer
T3	T3 Promoter Primer: 5'd[ATTAACCCTCACTAAAGGGA]3' -20mer
T7	T7 Universal Primer: 5'd[TAATACGACTCACTATAGGG]3' -20mer
T7T	T7 Terminator: 5'd[GCTAGTTATTGCTCAGCGG]3' -19mer